



Roots Vet Dental Supplies Ltd

Southfield Road Hinckley Leicestershire LE10 1UB, Tel:01455 698449

EMS/Woodpecker and DTE/Satelec Scaler tips

- G1 - Removal of supragingival deposits in all quadrants.
- G2 - Removal of heavy supragingival deposits on anterior teeth. Apply flat end to surface of teeth.
- G3 - Removal of supragingival deposits in all quadrants, including the interproximal and sulcus areas.
- G4 - Recommended for the treatment of interproximal spaces and for supragingival scaling.
- G5 - Recommended for treating simple cases and gross supragingival scaling.
- G6 - Recommended for removing voluminous supragingival deposits. Apply flat end to surface of teeth.
- P1 - Tips with thin points, recommended for cleaning and irrigation of subgingival deposits on root surfaces.
- P3 - Recommended for cleaning and irrigation of periodontal deep pockets.
- P4 - Recommended for cleaning and irrigation of periodontal shallow pockets.